MiniMed™ Quick-set™

Medtronic

Made in Mexico • Fabriqué au Mexique • Fabricado en México • Fabricado no México • 墨西哥制造 • メキシコ製

Medtronic MiniMed

18000 Devonshire Street
Northridge, CA 91325
USA
1 800 646 4633
+1 818 576 5555
www.medtronicdiabetes.com

Unomedical
A Convatec Company

Unomedical a/s
Aaholmvej 1-3, Østted
4320 Lejre
Denmark
www.infusion-set.com


(01)00763000400361

2020-06

M981717A6AA2_A

EN A. MiniMed tubing connector B. Tubing C. Needle protector feature D. Paper backing E. Site connector F. Adhesive tape G. Needle guard H. Disconnect cover (please save for later use)

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
Medtronic

EN Refer to the pump and reservoir user guides for installation and filling instructions.

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EN Fill the empty space in the soft cannula with the following volume of U-100 Insulin:

	6mm	9mm
FR	0.3u	0.5u
ES	0.3u	0.5u
PT	0.3e	0.5e
ZH	0.003ml	0.005ml
JA	0.003ml	0.005ml

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

Taiwan:
Medtronic (Taiwan) Ltd.
Tel: 02-21836000
Toll Free: +866-800-005285

Thailand:
Medtronic Thailand Ltd.
Tel: +662 232 7400

EN Disconnecting 1. Press the sides of the site connector. 2. Turn the site connector counterclockwise until the two arrows are aligned. 3. Lift up the site connector to remove it from the coupling site. 4. Connect the disconnect cover to the coupling site.

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

EN Reconnecting 1. Remove the disconnect cover from the coupling site, by pressing the grips, to reconnect the site connector. 2. Gently connect the site connector to the coupling site. 3. Turn the site connector clockwise until the arrow is aligned with the line on the adhesive tape.

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English

Indications for use

The MiniMed™ Quick-set™ is indicated for subcutaneous infusion of insulin from an infusion pump. The infusion set is indicated for single use.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

Contraindications

The infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous (IV) infusion. Do not use the infusion set for infusion of blood or blood products.

Warnings

Read all instructions carefully before using the infusion set. Failure to follow instructions may result in pain or injury.

Do not use the infusion set if the package is opened or damaged. An opened or damaged package may be contaminated. This can cause an infection. (For Medtronic pump users, do not use the infusion set if the tubing connector needle has been damaged.)

If using the infusion set for the first time, do so in the presence of a healthcare professional.

Inaccurate medication delivery, site irritation or infection may result from improper insertion or maintenance of the infusion site.

Do not reinsert the introducer needle into the infusion set. Reinsertion may cause tearing or damage of the soft cannula, which may result in unpredictable medication flow.

Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.

Remove the needle guard before inserting the infusion set.

Before insertion, clean the insertion site with a disinfectant.

If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

Precautions

Never prime or try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.

Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult pump manufacturer's instructions for use.

MiniMed™ Quick-set™ is a single-use item, which must be disposed of after usage. Do not clean or re-sterilize.

Check the infusion set frequently to ensure that the soft cannula remains firmly in place. Replace with a new infusion set at a new site, if the infusion set is not in place. Since the cannula is soft, it will not cause any pain if it slips out, and this may happen without notice. The soft cannula must always be completely inserted to receive the full amount of medication.

Do not put disinfectants, perfumes, or deodorants on the infusion set as these may impact the integrity of the infusion set.

Replace the infusion set every two to three days, or per healthcare professional's instructions.

If the infusion site becomes irritated or infected, replace the infusion set, and use a new site until the first site has healed. Replace the infusion set if the tape becomes loose, or if the soft cannula becomes fully or partially dislodged from the skin.

Use aseptic techniques when temporarily disconnecting the infusion set. Consult a healthcare professional on how to compensate for missed medication when disconnected.

Check blood glucose level one to three hours after insertion of the MiniMed™ Quick-set™ infusion set. Measure blood glucose on a regular basis. Discuss this with a healthcare professional.

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Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours after insertion.

If blood glucose level becomes inexplicably high, or an occlusion alarm occurs, check for clogs and/or leaks. If in doubt, change the infusion set, because the soft cannula may be dislodged or partially clogged. Discuss a plan with a healthcare professional to bring down blood glucose level. Test blood glucose level to make sure the problem is corrected.

Carefully monitor blood glucose level when disconnected and after reconnecting.

If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the national authority.

Storage and disposal

Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.

Store and handle insulin according to the manufacturer's instructions.

Dispose of the introducer needle in a proper sharps container according to local laws, to avoid the risk of needle stick.

Dispose of a used infusion set according to local regulations for biohazardous waste.

Warranty

For product warranty information, contact a local Medtronic support representative, or visit <http://www.medtronicdiabetes.com/warranty>.

This infusion set is compatible with the MiniMed™ Quick-setter™ insertion device. Refer to the insertion device user guide for insertion information.

© 2020 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (TM) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.

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